

Hospital ICU/CCU Electrical Bed

BA6000-N

- Material** : Mild steel
- Finishing** : Epoxy powder
- Safety load** : 230kg
- Backrest** : 0° - 70°
- Legrest** : 0° - 17°
- Kneerest** : 0° - 35°
- Trendelenburg & Reverse**
- Trendelenburg** : -15° to +15°
- Castor Dim(mm)** : 150mm dia.
- Mattress Dim(mm)** : 1980(L) x 915(W) x 100(H)
- Approx. Dim(mm)** : 2200(L) x 1050(W) x 480 - 780(H)

:: Product Features

- + HiLo, Backrest, Legrest, Trendelenburg & Reverse Trendelenburg position controlled by handset
- + PU side rail with split version
- + Caregiver control at foot end
- + X-ray translucent base at backrest
- + C/w CPR function
- + PE Head & Foot board
- + PU bumpers at four corners
- + IV drip pole storage and 6 IV pole attachment

:: Electrical System

- + Linak DC-actuator system (TUV & UL approval & EC60601-1)
- + 24V DC with built in backup battery
- + Voltage: 220 - 240V; 50Hz - 60Hz

:: Accessories Included

- + Hospital Medi-Flex Mattress
- + Telescopic IV Drip Pole with 2 hooks

:: Optional Accessories

- + Monkey pole
- + Oxygen cylinder holder
- + Portable monitor/Defibrillator tray
- + Orthopedic frame
- + Hospital Mattress with memory foam

Position:

