

ELECTRICAL BED
MODEL BA7001



Features

- Ultra low with the minimum height at 240mm (9.5” Inches)
- All parts are made by tooling to ensure uniform quality
- Finishing with epoxy powder coated & baked for chip and scratch resistance.
- Electric motorize bed with smooth and control operation by 3 actuators.
- Detachable head and foot boards
- Mobile on 4 x 100mm diameter castors with 2 brakes diagonally fitted
- Provide with drainage bag holder at both side of the bed
- Hilo, backrest & kneerest controlled by hand set
- Provision for IV pole
- With 4 corner rotating bumper

Technical Specifications

- Hilo: 240-650mm
- Working capacity: 180kg
- Backrest: 0° - 75°
- Auto Contour
- Kneerest: 0° - 45°
- Mattress Dim(mm): 1905(L) X 864(W) x 102(H)
- Overall Dim(mm): 2115(L) x 1013(W)

Electrical System

- Linak DC-actuator system (TUV & UL approved & EC60601-1)
- 24 DC Output
- Voltage: 220-240; 50Hz - 60Hz
- Built in backup battery

Optional

- Monkey/Lifting Pole (BX910)
- IV pole with 2 hooks- chrome plated (BX300)
- IV pole with 2 hooks- stainless steel (BX301)