


型号: KJW-D336(LN)

品名: 三功能电动护理床

THREE FUNCTIONS ELECTRIC CARE BED



主要功能:

- A:背部升降
- B:腿部升降
- C:整体升降

Main functions:

- A:Back section up/down.
- B:Knee section up/down.
- C:Whole bed up/down.

产品参数:

床体规格: L2130 × W980 × H430-690mm
 背部倾斜: $80^{\circ} \pm 5^{\circ}$
 腿部倾斜: $40^{\circ} \pm 5^{\circ}$
 床体可调高度: 430-690mm

Related parameters:

Overall dimension: L2130 × W980 × H430-690mm
 Back rest tilting: $80^{\circ} \pm 5^{\circ}$
 Knee board tilting: $40^{\circ} \pm 5^{\circ}$
 Hi-lo adjustment between: 430-690mm

标准配置:

高级医用电机三个, 控制盒一个;
 控制手柄一个, 豪华轮四个;
 塑料床头尾板一对, 弯管护栏一对;
 引流挂钩四个, 不锈钢输液杆一支;

Standard configuration:

3pcs Advanced medical actuators;
 1pc control box, 1pc handset;
 4pcs luxurious castors;
 1 pair plastic head/foot-end boards;
 1 pair L type siderails;
 4pcs urine hooks; 1pc S.S I.V pole;

可选配:

后配电源UPS; 快速放倒功能CPR;
 医用床垫; 餐板; 床头柜; 病历卡;
 杂物篓。

Options:

Backup power system (UPS);
 CPR function; Medical mattress;
 Dining board; Bedside cabinet;
 Patient file panel; Sundries basket.