



KAREWAY

广东康进威医疗科技有限公司

01

广东康进威医疗科技有限公司是一家专业致力于医疗器械的研发，生产，销售于一体的高新科技企业。本公司的产品以医用护理设备，康复设备为主，主要包括电动护理床，手动护理床，电动轮椅，手动轮椅，拐杖，助行器等。本公司自成立以来，汇聚了一批专业素质过硬，充满活力的工程技术人员。我们一直坚持“诚信，创新，服务”为理念，提供一流的产品，为客户创造价值。我们除了坚持独特的研发技术外，还特别坚持完善的产品售后服务体系。

本厂通过ISO9001:2008、ISO13485、CE 认证。我们致力于用创新和人性化的设计为客户提供完全，实用，便捷的产品，期待与您的合作。



Guangdong kareway medical technology Co., LTD is a professional medical equipment manufacturer with its own R&D, production and sales facility in Foshan China. Our company's products covers a range of electric hospital beds, manual hospital beds, electric wheelchairs, manual wheelchairs, walkers, walking sticks and other rehabilitation devices; we have a very professional R&D team with advanced manufacturing skills and good quality control system, after sale service is also our focus.

Our company has been certified by ISO9001:2008, ISO13485:2003. and our products have been approved by CE certificate. Our vision is to provide the world class products and services to attract the customers from all over the world. we look forward to good cooperation with you.

KAREWAY



KAREWAY

广东康进威医疗科技有限公司

KAREWAY



QUALIFICATION



TABLE OF CATALOGUE

01

ELECTRIC HOSPITAL BED
P07-P22

02

MANUAL HOSPITAL BED
P24-P36

03

OPERATING LIGHTS & TABLE SERIES
P38-P45

04

EMERGENCY STRETCHER
P46-P48

05

EMERGENCY TROLLEY
P50-P52

06

OTHER MEDICAL EQUIPMENT
P54-P59

07

WHEELCHAIR
P61-P65

08


REHABILITATION PRODUCTS
P66-P83

ELECTRIC HOSPITAL BED



KJW-D602

ELECTRIC BED WITH TWO COLUMN MOTORS



Main functions:

- A: Back section up/down.
- B: Knee section up/down.
- C: Whole bed up/down.
- D: Trendelenburg and Reverse Trendelenburg.

Related parameters:

Overall dimension : L2250xW1080xH520-910mm
 Back rest tilting: $80^{\circ} \pm 5^{\circ}$,
 Knee board tilting: $40^{\circ} \pm 5^{\circ}$
 Hi-lo adjustment between: 520-910mm;
 Reverse and trendelenburg $15^{\circ} \pm 2^{\circ}$;

Standard configuration:

2pcs column motors, 2pcs advanced medical actuators,
 1pcs control box,
 1pcs ACP (Attendant Control Operation Panel);
 central-controlled brake system,
 4pcs luxurious central locking castors;
 1 pair plastic head/foot-end board;
 1 pair Luxurious PP siderails
 4pcs urine hooks; 1pcs S.S I.V pole;


Options:

Backup power system (UPS);
 CPR function;
 Medical mattress; dining board;
 Bedside cabinet;



KJW-D501(PZ-D)

FIVE FUNCTIONS ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg;

Related parameters:

Overall dimension: L2250xW1080xH450-720mm;
 Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
 Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
 Hi-lo adjustment between: 450-720mm;
 Reverse and trendelenburg $15^{\circ} \pm 2^{\circ}$;

Standard configuration :

4pcs Advanced medical actuators,
 1pc control box, 1pc handset;
 Central-controlled brake system,
 4pcs luxurious central locking castors;
 1 pair plastic head/foot-end boards;
 1 pair Luxurious PP siderails
 4pcs urine hooks;
 1pc S.S I.V pole;


Options:

Backup power system (UPS);
 CPR function;
 Medical mattress;
 Dining board;
 Bedside cabinet;



KJW-D504(PZR)

FIVE FUNCTIONS ELECTRIC CARE BED
WITH BUILT IN CONTROL PANEL



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg;

Related parameters:

- Overall dimension: L2300xW1050xH480-790mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 480-790mm;
- Reverse and trendelengburg $15^{\circ} \pm 2^{\circ}$;

Standard configuration :

- 4pcs advanced medical actuators;
- 1pc control box ;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end board with built in control panel ;
- 1 pair PP siderails with built in control panels;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D309(PZR)

THREE FUNCTIONS ELECTRIC CARE BED
WITH BUILT IN CONTROL PANEL



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:


- Overall dimension: L2200xW1060xH480-760mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 480-760mm;

Standard configuration :

- 3pcs advanced medical actuators;
- 1pc control box ;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards with built in control panel;
- 1 pair PP siderails with built in control panels;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D307(PZ-S1)

THREE FUNCTIONS ELECTRIC&MANUAL CARE BED



Main functions:

- A:Back section up/down;
- B:Knee section up/down;
- C:Whole bed up/down;
- D:Electric/manual operation;

Related parameters:

- Overall dimension: L2170xW990xH460-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 460-720mm;

Standard configuration:

- 3pcs advanced medical actuators;
- 1pc control box; 1pc handset;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair Luxurious PP siderails;
- 4pcs urine hooks;1pc S.S I.V pole;



Options:

- Backup power system (UPS);
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D307(PZ-S2)

THREE-FUNCTION ELECTRIC MANUAL NURSING BED (LENGTHENED AND WIDENED)



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

Overall dimension: L2230 × W1140 × H480–760mm;
 Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
 Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
 Hi-lo adjustment between: 480–760mm;

Standard configuration:

3pcs advanced manual&electric medical actuators;
 1pc control box, 1pc handset;
 central-controlled brake system;
 4pcs luxurious central locking castors;
 1 pair plastic head/foot-end boards;
 1 pair plastic siderails;
 4pcs urine hooks;1pc S.S I.V pole;1pc handset;


Options:

Backup power system (UPS);
 Medical mattress; Dining board; Bedside cabinet;



KJW-D504(PZ-C)

FIVE FUNCTIONS ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg;

Related parameters:


Overall dimension: L2300xW1050xH480-790mm;
 Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
 Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
 Hi-lo adjustment between: 480-790mm;
 Reverse and trendelenburg $15^{\circ} \pm 2^{\circ}$;

Standard configuration:

4pcs advanced medical actuators;
 1pc control box; 1pc handset;
 Central-controlled brake system,
 4pcs luxurious central locking castors;
 1 pair plastic head/foot-end boards;
 1 pair PP siderails;
 4pcs urine hooks; 1pc S.S. I.V pole;


Options:

Backup power system (UPS);
 CPR function;
 Medical mattress;
 Dining board;
 Bedside cabinet;



KJW-D502(PN-C)

FIVE FUNCTIONS ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg;

Related parameters:


- Overall dimension: L2170 × W990 × H430-690mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430-690mm;
- Reverse and trendelenburg $15^{\circ} \pm 2^{\circ}$;

Standard configuration :

- 4pcs advanced medical actuators;
- 1pc control box; 1pc handset;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair PP siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D303(PZ-C)

THREE FUNCTIONS ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170xW990xH460-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 460-720mm;

Standard configuration :

- 3pcs advanced medical actuators;
- 1pc control box; 1pc handset;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair PP siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D302(LZ-C)

THREE FUNCTIONS ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170xW990xH430-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430-720mm;

Standard configuration:

- 3pcs advanced medical actuators;
- 1pc control box, 1pc handset;
- Central-controlled brake system;
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D301(LN-C)

THREE FUNCTIONS ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170xW990xH430-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430-720mm;

Standard configuration :

- 3pcs advanced medical actuators;
- 1pc control box, 1pc handset;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D305(LN-C)

THREE FUNCTIONS ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170xW990xH290-680mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 290-680mm;

Standard configuration :

- 3pcs advanced medical actuators;
- 1pc control box, 1pc handset;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D201(LN-C)

TWO FUNCTIONS ELECTRIC CARE BED



Main functions:

- A:Back section up/down;
- B:Knee section up/down;

Related parameters:

- Overall dimension : L2170xW990xH490mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;

Standard configuration :

- 2pcs advanced medical actuators;
- 1pc control box,1pc handset;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D503(LN)

TURNING OVER ELECTRIC CARE



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Electric turning over bed;

Related parameters:

Overall dimension: L2120xW940xH560mm;

Standard configuration:

- 4pcs advanced medical actuators;
- 1pc control box; 1pc handset;
- 4pcs luxurious double-face castors;
- 1 pair PP head/foot-end boards;
- 1 pair aluminium siderails;
- 4pcs urine hooks;
- 1pc S.S. I.V pole;



Options:

- Backup power system (UPS);
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D331(LN)

THREE FUNCTIONS ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2100xW970xH430-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430-720mm;

Standard configuration:

- 3pcs advanced medical actuators;
- 1pc control box;
- 1pc handset;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



Options:

- Backup power system (UPS);
- CPR function;
- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-D335(LN)

THREE FUNCTIONS SUPER LOW ELECTRIC CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

Overall dimension: L2100 × W970 × H270–680mm;
 Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
 Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
 Hi-lo adjustment between: 270–680mm;

Standard configuration:

3pcs advanced medical actuators;
 1pc control box;
 1pc handset;
 4pcs luxurious castors;
 1 pair plastic head/foot-end boards;
 1 pair aluminium siderails;
 4pcs urine hooks;
 1pc S.S I.V pole;

Options:

Backup power system (UPS);
 CPR function;
 Medical mattress;
 Dining board;
 Bedside cabinet;

MANUAL HOSPITAL BED

A photograph of a family celebrating a patient's recovery in a hospital room. A young girl with blonde hair is sitting up in a hospital bed, wearing a light blue hospital gown. She is smiling and looking towards a young boy sitting on the edge of the bed. The boy is wearing a striped shirt and is also smiling. A woman in a pink shirt and a man in a purple polo shirt are sitting on the bed with them, all looking at a card or document. A teddy bear is sitting on the bed next to the girl. There are balloons and flowers in the background, suggesting a birthday or a special occasion. The scene is set in a bright, clean hospital room with a window in the background.



KJW-S403(LZ)

MANUAL HYDRAALIC HOSPITAL BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg

Related parameters:

- Overall dimension: L2170xW1080xH460-750mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 460-750mm;
- Trendelenburg: $12^{\circ} \pm 2^{\circ}$;
- Reverse trendelenburg: $6^{\circ} \pm 2^{\circ}$

Standard configuration :

- Manual Hydraulic systems;
- Central controlled brake system;
- 4pcs luxurious central locking castors;
- 1pair plastic head/foot-end boards;
- 1pair L type sidrails;
- 4pcs urine hooks;
- 1pcs S.S LV pole


Options:

- Medical mattress;
- Dining boird;
- Bedside cabiner;



KJW-S332(SN)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: With food table

Related parameters:

- Overall dimension: L2170xW980xH450-885mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 450-885mm;

Standard configuration :

- 3 sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 4 pcs S. S side rails. ;
- 4pcs urine hooks;
- 1pc S.S I.V pole;

Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-S308(LN-B)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:


- Overall dimension: L2160xW920xH480-760mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 480-760mm;

Standard configuration:

- 3 set manual crank systems;
- 4 pcs luxurious castors, 4 pcs legs without castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4 pcs urine hooks;
- 1 pcs S.S.IV pole;


Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-S402(PZ)

FOUR CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg;

Related parameters:

- Overall dimension: L2170 × W990 × H460–750mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 460–750mm;
- Reverse and trendelenburg $15^{\circ} \pm 2^{\circ}$;

Standard configuration :

- 4 sets manual crank systems;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-S401(LN-C)

FOUR CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg;

Related parameters:

- Overall dimension: L2170 × W990 × H430–720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430–720mm;
- Reverse and trendelenburg $15^{\circ} \pm 2^{\circ}$;

Standard configuration:

- 4 sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-S303(PZ-C)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170 × W990 × H430–720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430–720mm;

Standard configuration:

- 3 sets manual crank systems;
- Central-controlled brake system;
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair Luxurious PP siderails;
- 4pcs urine hooks;
- 1pc S.S. I.V pole;


Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;
- Sundries basket;



KJW-S301(LN-C)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170xW990xH430-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430-720mm;

Standard configuration :

- 3 sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;
- Sundries basket;



KJW-S302(LZ-C)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170xW990xH430-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430-720mm;

Standard configuration:

- 3 sets manual crank systems;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;
- Sundries basket;



KJW-S202(PN-C)

TWO CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170 × W990 × H480mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;

Standard configuration :

- 2 sets manual crank systems;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;
- Sundries basket;



KJW-S201(LN-C)

TWO CRANKS MANUAL CARE BED



Main functions:

A:Back section up/down;
B:Knee section up/down;

Related parameters:

Overall dimension : L2170xW990xH480mm;
Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
Knee board tilting: $40^{\circ} \pm 5^{\circ}$;

Standard configuration :

2 sets manual crank systems ;
4pcs luxurious castors ;
1 pair plastic head/foot-end boards ;
1 pair L type siderails ;
4pcs urine hooks ;
1pc S.S I.V pole ;



Options:

Medical mattress ;
Dining board ;
Bedside cabinet ;
Sundries basket ;



KJW-S131(L)

SINGLE CRANK MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;

Related parameters:

Overall dimension:
L2130 x W980 x H490mm;

Standard configuration:

- 1 set manual crank systems;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails;



Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



KJW-S131(LN)

SINGLE CRANK MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;

Related parameters:


Overall dimension:
L2100xW970xH480mm;
Back rest tilting: $80^{\circ} \pm 5^{\circ}$;

Standard configuration:

- 1 set manual crank systems;
- 4 pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails;
- 4 pcs urine hooks;
- 1 pc S.S I.V pole;

Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



KJW-S231(LN)

TWO CRANKS MANUAL CARE BED

Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:



- Overall dimension:
L2100 × W970 × H480mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;

Standard configuration :

- 2 sets manual crank systems ;
- 4pcs luxurious castors ;
- 1 pair plastic head/foot-end boards ;
- 1 pair Aluminium siderails
- 4pcs urine hooks ;
- 1pc S.S I.V pole ;

Options:

- Medical mattress ;
- Dining board ;
- Bedside cabinet ;
- Sundries basket ;



KJW-S331(LN)

THREE CRANKS MANUAL CARE BED

Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension:
L2100 × W970 × H360-680mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between:
360-680mm;

Standard configuration :

- 3 sets manual crank systems ;
- 4pcs luxurious castors ;
- 1 pair plastic head/foot-end boards ;
- 1 pair aluminium siderails
- 4pcs urine hooks ;
- 1pc S.S I.V pole ;


Options:

- Medical mattress ;
- Dining board ; Bedside cabinet ;



KJW-S331(G)

THREE CRANKS MANUAL ORTHOPEDICS BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;

Related parameters:

- Overall dimension :
- L2220xW920xH480mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;

Standard configuration:

- 3sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails
- 4pcs urine hooks;
- 1pc S.S I.V pole;
- 1set orthopedics traction frame;


Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



KJW-S432(G)

FOUR CRANKS MANUAL ORTHOPEDICS BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension :
- L2220xW920xH460-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between:
- 460-720mm;

Standard configuration:

- 4sets manual crank systems;
- control-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails
- 4pcs urine hooks;
- 1pc S.S I.V pole;
- 1set orthopedics traction frame;

Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;

OPERATING LIGHTS & TABLE SERIES

Shadowless lamp series

KJW-OPL20/10



KJW-OPL20



KJW-OPL10



PWM dimming technology,
no flickering.

LED cold light,
head temperature rise
less than 1°C



Germany brand Ondal spring arms



Removable sterilizer handle
Light field adjusting



OSRAM LED



Beautiful ultra-thin
lamp

Type	OPL20 Main lamp	OPL10 Main lamp
Head Diameter	700mm	500mm
LEDS		
Illumination		
CRI		
Color Temperature		
Illumination Depth		
Spot Diameter		
LED Life		
LED Brand		
UPS		

KJW-OPL10B

OPL10B standing operating light uses LED cold light source, it can be used for surgery and examination.



Type	PLM-700 Main lamp	PLM-500 Main lamp
Head Diameter	700mm	500mm
LEDS	80	48
Illumination	40~160klx	
CRI	Ra > 95	
Color Temperature	3500~5000K	
Illumination Depth	1200mm	
Spot Diameter	120~300mm	
LED Life	> 50,000Hours	
LED Brand	OSRAM	
Battery	Standard	



LCD control panel

10 steps adjustable illumination

5 steps adjustable color temperature



Germany brand Ondal spring arms



Big capacity battery Castors with brakes

KJW-D9600

Electric Hydraulic Operating table



D9600 is electric hydraulic operating table, which all functions are electric hydraulic controlled.



Four point electric hydraulic brake
 One button zero position
 Leg section is also electric hydraulic controlled

Sliding distance 350cm, which is C-arm compatible
 Equipped with side controller

With flex/reflex function
 Big capacity battery for 50-80 times surgery



Main specification

Length/Width	≥170mm/550mm	Back Section Up/Down	≥90°/≥40°
Height	690mm-1100mm	Sliding Distance	≥350mm
Trendelenburg/Revers Trendelenburg	≥30°/≥30°	Kidney Bridge Lift	110mm
Lateral Tilt	≥20°/≥20°	Flex/Reflex	≥90°/≥220°
Head Section up/down	≥40°/≥90°	Load Capacity	250kg
Head Section up/down/outside	≥15°/≥90°/≥90°		



KJW-D9500

Electric Operating table

D9500 is electric operating table, which all functions are controlled by imported motor.

Sliding distance is 350mm, which is C-arm compatible


8cm memory sponge mattress

Big capacity battery for 50-80 times surgery.



Main specification

Length/Width	≥2100mm/550mm	Leg Section Up/Down/outside	≥15°/≥90°/≥90°
Height	700mm-1000mm	Back section Up/Down	≥75°/≥20°
Trendelenburg/Revers Trendelenburg	≥25°/≥25°	Sliding Distance	≥350mm
Lateral Tilt	≥15°/≥15°	Kidney Bridge Lift	110mm
Head Section up/down	≥40°/≥90°	Load Capacity	250kg



Electric OR table series

KJW-OP-D9200A/B

Electric Operating table

OP-D9200A/B is electric operating table, which all functions are controlled by high quality medical motor

- Memory sponge mattress
- Big capacity battery for 50-80 times surgery



KJW-OP-D9200A

Sliding distance is 350mm, which is C-arm compatible



KJW-OP-D9200B

Main specification

Model	2000mm/500mm	2020mm/500mm
Length/Width	2000mm/500mm	2020mm/500mm
Height	750mm-1000mm	550mm-800mm
Trendelenburg/Revers Trendelenburg	≥25°/≥25°	≥25°/≥20°
Lateral Tilt	≥20°/≥20°	≥20°/≥20°
Head Section Up/Down	≥45°/≥90°	≥30°/≥90°
Leg Section Up/Down/Outside	≥15°/≥90°/≥90°	≥15°/≥90°/≥90°
Back Section Up/Down	≥70°/≥35°	≥45°/≥20°
Sliding Distance	≥350mm	—
Kinney Bridge Lift	100mm	80mm
Load Capacity	250kg	250kg



Skull clamp compatible



KJW-OP-S8200
Hydraulic Operating table

KJW-OP-S8200 is hydraulic mechanical operating table. The table up-down is hydraulic controlled



KJW-OP-S8300
Hydraulic Operating table

KJW-OP-S8300 is hydraulic mechanical operating table. The table up-down is hydraulic controlled



The column and base cover are all stainless steel 304#, which is convenient for clean and disinfection.

Main specification

Length/Width	2080mm/480mm	Leg Section Up/Down/Outside	≥15°/≥90°≥90°
Height	750mm-1000mm	Back Section Up/Down	≥90°/≥10°
Trendelenburg/Revers Trendelenburg	≥25°/≥25°	Kidney Bridge Lift	120mm
Lateral Tilt	≥18°/≥18°	Load Capacity	250kg
Head Section up/down	≥45°/≥90°		

Main specification

Length/Width	2000mm/480mm	Leg Section Up/Down/Outside	≥15°/≥90°≥90°
Height	720mm-980mm	Back Section Up/Down	≥75°/≥18°
Trendelenburg/Revers Trendelenburg	≥30°/≥20°	Kidney Bridge Lift	115mm
Lateral Tilt	≥18°/≥18°	Load Capacity	250kg
Head Section up/down	≥45°/≥90°		



The table top are all stainless steel 304#, which is durable and convenient for clean and disinfection



Electric obstetric delivery bed series

KJW-OP-D9400
Electric Obstetric table

KJW-OP-D9400 electric delivery bed up-down, back up-down, seat section up-down (or Trendelenburg and rever-trendelenburg) are all controlled by

- ▮ The mattress is seamless foaming pad, which is beautiful and easily for clean.
- ▮ The leg board can be opened 90 ° , which is convenient for operation.
- ▮ The foot controller is standard configuration.



Main specification

Length/Width	1900mm/600mm	Back Section Up/Down	≥70°/≥10°
Height	600mm-900mm	Leg section opened	≥90°/≥10°
Trendelenburg/Revers Trendelenburg	≥20°/≥20°	Load Capacity	250kg



All the cover, able top, basin are all made of stainless steel 304#, which is convenient for clean and disinfection.

The examination light is optional, which can be used for gynecology examination.

Hydraulic obstetric delivery bed

KJW-OP-S8400
Hydraulic Obstetric table

KJW-OP-S8400 delivery bed table up-down is hydraulic controlled.

The column and base cover, basin are all made of stainless steel 304#, which is convenient for clean and disinfection.



Main specification

Length/Width	2000mm/600mm
Height	620mm/860mm
Trendelenburg/Revers Trendelenburg	≥20°/≥20°
Back Section Up/Down	≥65°/≥10°
Leg Section Down	250kg

Manual obstetric delivery bed

KJW-OP-S8200
Hydraulic Obstetric table



Main specification

Length/Width/Height	1900mm/600mm/750mm
Back Section Up/Down	≥45°/≥10°
Seat Section Up/Down	≥25°/≥8°
Leg Section Up/Down	≥25°/≥25°

KJW-OP-S8200 Ordinary obstetric bed is used for obstetrics gynecology examination. The bed top and main structure are all made of stainless steel 304#

EMERGENCY STRETCHER

KJW-ES003

CONNECTING STRETCHER

SIZE: L3650×W640×H640-970mm

Standard Devices:

- Two manual cranks;
- one back rise/fall gas spring;
- Lockable and tuckaway side rail;
- Luxurious central locking castors;
- 5th wheel;
- Oxygen cylinder holder;
- one IV pole;
- Two mattresses;





KJW-ES001

LUXURIOUS MANUAL EMERGENCY STRETCHER

SIZE: L1930×W640×H450-840mm

Standard Devices:

- One manual crank;
- one back rise/fall gas spring;
- Lockable and tuckaway side rail;
- Four luxurious central locking castors;
- 5th wheel;
- one oxygen cylinder holder;
- one IV pole;
- one mattress;



KJW-ES005

LUXURIOUS HYDRAULIC EMERGENCY STRETCHER

SIZE:L2080 x W850 x H700~1000mm

Standard Devices:

- One manual crank;
- Two back rise/fall gas springs;
- Fold down side rail;
- Two hydraulic actuator;
- Four central locking castors;
- 5th wheel;
- one oxygen cylinder holder;
- one IV pole;
- one mattress;



KJW-ES006

LUXURIOUS HYDRAULIC EMERGENCY STRETCHER

SIZE: L1930*W750*H700-980mm

Standard Devices:

- Two back rise/fall gas springs;
- Lockable and tuckaway side rail;
- Two hydraulic pumps;
- Four central locking castors;
- 5th wheel;
- one oxygen cylinder holder;
- one IV pole;



KJW-ES008

LUXURIOUS HYDRAULIC EMERGENCY STRETCHER

SIZE:L2160*W780*H560-890mm

Standard Devices:

- back rise/fall gas spring;
- Lockable and tuckaway side rail;
- Two hydraulic pumps;
- Four central locking castors;
- 5th wheel;
- one oxygen cylinder holder;
- one IV pole;
- one mattress;
- One writing table;

EMERGENCY TROLLEY





KJW-ET001

- Emergency Trolley
- SIZE: L:850 × 520 × 960mm
S:750 × 480 × 920mm



KJW-ET002

- Emergency Trolley
- SIZE: L:850 × 520 × 960mm
S:750 × 480 × 920mm



KJW-ET005

- Emergency Trolley
- SIZE: 650 × 480 × 920mm



KJW-ET006

- Emergency Trolley
- SIZE: 650 × 480 × 920mm



KJW-AT003

- Anesthesia Trolley
- SIZE: 750 × 480 × 920mm



KJW-AT004

- Anesthesia Trolley
- SIZE: 750 × 480 × 920mm



KJW-AT005

- Anesthesia Trolley
- SIZE: 650 × 480 × 920mm



KJW-AT006

- Anesthesia Trolley
- SIZE: 650 × 480 × 920mm



KJW-TT001

- Treatment Trolley
- SIZE: 750 × 480 × 920mm



KJW-TT002

- Treatment Trolley
- SIZE: 750 × 480 × 920mm



KJW-TT003

- Treatment Trolley
- SIZE: 750 × 480 × 920mm



KJW-TT004

- Treatment Trolley
- SIZE: 600 × 470 × 1090mm



KJW-MT001

- Medicine Trolley
- SIZE: 600 × 470 × 1090mm



KJW-MT012B

- Medicine Trolley
- SIZE: 750 × 480 × 920mm



KJW-MT013B

- Medicine Trolley
- SIZE: 750 × 480 × 920mm



KJW-MT014B

- Medicine Trolley
- SIZE: 750 × 480 × 920mm

**KJW-IT001**

- Infusion Trolley
- SIZE: 750 × 480 × 920mm

**KJW-IT002**

- Infusion Trolley
- SIZE: 750 × 480 × 920mm

**KJW-IT003**

- Infusion Trolley
- SIZE: 750 × 480 × 920mm

**KJW-MNT001**

- Morning Nursing Trolley
- SIZE: 1000 × 540 × 950mm

**KJW-MNT002**

- Morning Nursing Trolley
- SIZE: 1000 × 540 × 950mm

**KJW-LNT01T**

- Linen Trolley
- SIZE: 600 × 420 × 900mm

**KJW-LNT01S**

- Linen Trolley
- SIZE: 600 × 420 × 900mm

**KJW-LNT02T**

- Linen Trolley
- SIZE: 800 × 590 × 900mm

**KJW-LNT02S**

- Linen Trolley
- SIZE: 800 × 590 × 900mm



OTHER MEDICAL EQUIPMENT




KJW-M01

- Mattress
- Size: 1920 × 840 × 80mm
- Material:wearable cloth cover,sponge- cocount&plam fiber



KJW-M02

- Mattress
- Size: 1920 × 840 × 100mm
- Material:wearable PVC cover,high density sponge



KJW-M04

- Mattress
- Size: 1920 × 840 × 80mm
- Material:wearable PU cover,natural latex



KJW-BC604

- Plastic Baby Cart
- Size: 880 × 510 × 820-980mm



KJW-BC600

- Plastic Baby Cart
- Size: 810 × 510 × 870mm
- Material:stainless steel



KJW-BC601

- Plastic Baby Cart
- Size: 810 × 460 × 930mm
- Material:stainless steel

2.8" / 5.6" Bubble Pad Stage I Overlay System

KJW-AM02

Alternating Bubble Mattress With Pump System



- Recycled Parts
- Flame Retardant
- Water Resistant
- Odor Free

M7-303 is designed for the treatment and prevention of pressure sore-stage I in short-term and domiciliary therapies. The mattress in PVC is composed of 130 7cm-high bubble cells, particularly comfortable. It is fixed to the bed by extra flaps on both top and bottom side. Very easy to install, maintain and use thanks to the

- ✓ Without Flaps Optional
- ✓ Low Air Loss Holes Optional

Color Box Code: B2016Y1 (2.8")
20' Container: 5100 sets (Bubble Pad & Pump)
7100

possibility of adjust comfort range. The pump can be hung to 140" HQ: sets (Bubble Pad & Pump)



5.6" Bubble Mattress

Technical Parameters Of The Bubble Mattress (2.8")

Dimensions without Flaps	M7-303A: 200*90*7cm (Standard Size) M7-303B: 90*60*7cm (Economic Size)
Extension Flaps	Head & Feet Side
Mattress material	Softness Medical-Grade PVC
Material thickness	0.32mm
Cold-resistant	-30°C
Number of cells	130
Weight support	155kg
HF Welding	One-Time Forming Paper TOSHIBA Power Tube
Certificate	RoHS / CE-ET / EN 21
Warranty	2 months

Economic 4" Cell Stage II Overlay System

KJW-AM04

9.5cm / 4" Cell Alternating Pressure Mattress With Pump



Air Tube
Nylon+PVC

1 Five Spare Cell

- Recycled Parts
- Flame Retardant
- Water Resistant
- Odor Free

M6-601 Alternating Pressure Mattress is made up by 20 pcs of 9.5cm / 4" cells. The Alternating System includes A/B groups to keep inflating and deflating alternatively in cycles which forming the wave imitation. The alternating function of interchangeable cells with pumps can reduce the pressure part between the patient's body and mattress in order to prevent and relief the pressure ulcers.



Two tubes alternating



The 20 9.5cm-high cells are fixed by means of pins+staps. More easy to install on 3.5" connecting tubes.



Medical-Grade Nylon+PVC Thickness: 0.24mm

Air Pump (P9000)

the Most Quiet Pump in China



Technical Parameters Of The Mattress	9.5cm / 4" Cell
Dimensions (L*W*H)	220*100*100mm
Cell Height	9.5cm (4")
Cell and tube material	Nylon+PVC
System Color	Blue/Black
Number of Cells	20*4 Pins+staples
Test Duration	24 hours (Standard)
Pressure Loss Stage	1.5"
Mattress Material	0.32mm
Weight Support	155kg
Warranty	2 months

www.kareway.com support@kareway.com



KJW-BC02

- Abs Beside Cabinet
- Size: 480 × 480 × 760mm



KJW-BC04

- Bedside Cabinet
- Size: 490 × 460 × 930mm



KJW-BC06

- Stainless Steel Bedside Cabinet
- Size: 495 × 420 × 790mm



KJW-FT01

- Food Table



KJW-MT01

- Mobile Table
- Size: 820 × 450 × (700–900)mm



KJW-MT05

- Mobile Table
- Size: 900 × 450 × (800–1090)mm



KJW-SS100B

- Emergency Trolley
- Size: 1900*640*710mm
- Material: stainless steel
- Description: back adjustable, height adjustable



KJW-GEB02

- Gynecological Examination Couch
- Size: 1940 x 600 x 790mm
- Material: stainless steel with foam mattress
- Description: back adjustable, height adjustable



KJW-SS100C

- Emergency Trolley
- Size: 1800*620*(530-760)mm
- Material: stainless steel
- Description: back adjustable, height adjustable



KJW-GEB03

- Gynecological Examination Couch
- Size: 1540*600*790mm
- Material: stainless steel with foam mattress
- Description: back adjustable, height adjustable



KJW-WS04

- Stainless Steel Ward Screen
- Size:2400 × 1800mm



KJW-ET802

- Examination Table
- Size:1920*600*730mm
- Material:stainless steel
- Description:weight capacity 300KGS,back adjustable
- fixed hight



KJW-ET803

- Examination Table
- Size:1920*600*700mm
- Material:stainless steel
- Description:back adjustable, fixed hight
- With draws and cabinets with one step



KJW-ET801

- Examination Table
- Size:1900*600*700mm
- Material:stainless steel
- Description:back adjustable, fixed hight



KJW-ST02A

- Powder Coated Stool
- SIZE: 450 × 430 × 350mm



KJW-AC01

- Accompanying Chair
- Size: folded 650 × 610 × 950mm
unfolded 650 × 610 × 1900mm



KJW-IVS03

- I.V Stand
- Size: base: Φ 360mm
- I.V pole: 1200–2000mm
- Material: painted
- Metal steel base with plastic cover S.S pole



KJW-IVS07

- I.V Stand
- Size: base: Φ 520mm
- I.V pole: 1200–2000mm
- Material: S.S base
- S.S pole



KJW-IVS08

- I.V Stand
- Size: base: Φ 360mm
- I.V pole: 1200–2000mm
- Material: painted
- Metal steel base S.S pole



KJW-OC01

- Oxygen Cylinder Cart
- Size: 500 × 450 × 1200 mm



KJW-OC02

- Oxygen Cylinder Cart
- Size: 500 × 450 × 1200mm

WHEELCHAIR



KJW-813LA

ELECTRIC WHEELCHAIR

Frame Material: Aluminium;
 Armrest & Footrest: Flip up armrest ,detachable footrest;
 Seat width: 18"/20" inches;
 Front castors: 8" big PU tyre ,
 Rear wheels: 12" big PU tyre;
 Speed: 1~7km/h (adjustable);
 Battery: Lead-acid battery(2pcs) ;
 Per charge range: 18km;



KJW-811A

ELECTRIC WHEELCHAIR

Frame Material: Metal Steel;
 Armrest & Footrest: Flip up armrest ,detachable footrest;
 Seat width: 18" inches;
 Front castors: 8" big PU tyre ,
 Rear wheels: 12" big PU tyre;
 Speed: 1~7km/h (adjustable);
 Battery: Lead-acid battery(2pcs) ;
 Per charge range: 18km;



KJW-813LB

ELECTRIC WHEELCHAIR

Frame Material: Aluminium;
 Armrest & Footrest: Flip up armrest ,detachable footrest;
 Seat width: 18"/20" inches;
 Front castors: 8" big PU tyre ,
 Rear wheels: 24" big PU tyre;
 Speed: 1~7km/h (adjustable);
 Battery: Lead-acid battery(2pcs) ;
 Per charge range: 18km;



KJW-820B

ELECTRIC WHEELCHAIR

Frame Material: Metal Steel ;
 Armrest & Footrest: Detachable armrest, detachable footrest;
 Seat width: 17" inches;
 Front castors: 8" big PU tyre ,
 Rear wheels: 12" big PU tyre;
 Speed: 1~7km/h (adjustable);
 Battery: Lead-acid battery(2pcs) ;
 Per charge range: 18km;

型号	尺寸 (mm)											座高	座宽	座深
	总宽	扶手宽	扶手高	扶手距	扶手距	扶手距	扶手距	扶手距	扶手距	扶手距	扶手距	座高	座宽	座深
KJW-813LA	62	36	46	32	20	52	93	110	43	43	110	85X42X78	45	50
KJW-813LB	67	42	46	61	20	52	96	110	43	43	110	86X43X80	45	50

型号	尺寸 (mm)											座高	座宽	座深
	总宽	扶手宽	扶手高	扶手距	扶手距	扶手距	扶手距	扶手距	扶手距	扶手距	扶手距	座高	座宽	座深
KJW-811A	46	33	46	32	20	51	94	103	40	46	110	80X38X76	56	61
KJW-820B	64	-	46	32	20	53	107	127	43	74	110	76X44X79	65	70



KJW-616(LZJ-L)

- 1) Aluminum lightweight frame
- 2) Fixed armrest and footrest
- 3) 6" solid castor and 20"PU rear wheel with aluminum frame
- 4) Drop back handle with united brake



KJW-666(LZJ-S)

- 1) Aluminum lightweight frame
- 2) Drop back handle with unitell brake
- 3) Flip up armrest and footrest
- 4) Solid 6" pvc castor and 16"pu back wheel



KJW-613(LZJ-L)

- 1) Aluminum lightweight frame; double seat and back (net seat + detachable soft cushion) inner net seat for summer use, outer detachable cushion for winter use
- 2) Drop back handle with united brake, fixed armrest and footrest
- 3) 8" PVC castor, 24" PU rear wheel with aluminum frame



KJW-616(LZJ-S)

- 1) Aluminum lightweight frame
- 2) Drop back handle with united brake, fixed armrest and footrest
- 3) 6" solid castor and 12" PU rear wheel



KJW-666(LZJ-L)

- 1) Aluminum lightweight frame
- 2) Drop back handle with united brake
- 3) Flip up armrest and detachable footrest
- 4) Solid 6" PVC castor and 24"PU rear wheel



KJW-613(LZJ-S)

- 1) Aluminum lightweight frame
- 2) Fixed armrest and footrest
- 3) 8" solid castor and 24" PU mag rear wheel

型号	尺寸 (cm)	座宽	座深	座高	脚踏高	脚踏宽	脚踏深	脚踏距	脚踏高	脚踏宽	脚踏深	脚踏距	脚踏高	脚踏宽	脚踏深	脚踏距	脚踏高	脚踏宽	脚踏深	脚踏距	脚踏高	脚踏宽	脚踏深	脚踏距	
KJW-616LZJ-L	57	25	41	50	15	49	88	86	36	39	100	76X25X71	11	13											
KJW-666LZJ-S	65	29	46	40	15	47	89	98	40	40	120	72X29X72	11	13											
KJW-613LZJ-L	62	24	46	59	20	50	88	97	40	42	100	97X24X76	11	13											

型号	尺寸 (cm)	座宽	座深	座高	脚踏高	脚踏宽	脚踏深	脚踏距	脚踏高	脚踏宽	脚踏深	脚踏距	脚踏高	脚踏宽	脚踏深	脚踏距	脚踏高	脚踏宽	脚踏深	脚踏距	脚踏高	脚踏宽	脚踏深	脚踏距	
KJW-616LZJ-S	52	22	41	30	15	49	89	85	36	39	75	68X22X71	11	13											
KJW-666LZJ-L	65	30	46	50	15	47	89	102	40	40	120	77X30X72	12	14											
KJW-613LZJ-S	62	24	46	59	20	50	86	97	40	42	100	97X24X86	11	13											



KJW-615(LZJ-S)

- 1) Collapsible portable aluminum frame, with united brake
- 2) Drop back handle, fixed armrest and tip-back footrest
- 3) 6" solid castor and 12" PU rear wheel



KJW-615(LZJ-L)

- 1) Aluminium frame,withunited brake
- 2)Drop back handle,fixed armrest and tip-back footrest
- 3)6"solid castor and 20"PU rear wheel with aluminium frame



KJW-662(LZJ-S)

- 1) Aluminum frame,with united brake
- 2) Drop back handle, flip up armrest and detachable footrest
- 3) 6" solid castor and 16" PU back wheel.



KJW-662(LZJ-L)

- 1) Aluminum frame,with united brake
- 2) Drop back handle, flip up armrest and detachable footrest
- 3) 6" solid castor and 24" PU back wheel.



KJW-635L

- 1) Aluminum frame
- 2) Flip up armrest and detachable footrest
- 3) 8" solid castor and 24" PU back wheel.



KJW-715LHE

- 1)Aluminum frame
- 2)Detachable armrest and footrest
- 3) 6" solid castor and 16" PU rear wheel.

型号	尺寸 (mm)											U1	U2	U3	U4	U5	U6	U7	U8	U9	U10	U11	U12		
KJW-615LZJ-S	64	29	46	30	15	49	98	102	42	45	100	71X31X70	11	13											
KJW-662LZJ-S	64	29	46	40	15	49	98	102	42	45	100	71X31X70	12	14											
KJW-635L	64	29	46	40	15	49	98	102	42	45	100	71X31X70	12	14											

型号	尺寸 (mm)											U1	U2	U3	U4	U5	U6	U7	U8	U9	U10	U11	U12		
KJW-615LZJ-L	62	27	46	50	15	46	88	100	36	42	100	77X33X72	11	13											
KJW-662LZJ-L	62	27	46	60	15	46	88	100	36	42	100	77X33X72	13	15											
KJW-715LHE	58	33	46	40	15	43	96	112	40	45	100	86X33X100	29	32											



KJW-611

- 1) Chromed steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor, 24" solid rear



KJW-636M

- 1) Powder coated metal steel frame
- 2) Flip up armrest and detachable footrest
- 3) 8" solid castor and 24" solid back



KJW-621

- 1) Chromed steel frame
- 2) Fixed armrest and detachabl footrest
- 3) 8" PVC front castor, 24" solid rear wheel



KJW-635(M)

- 1) Metal steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 200X50 PVC front castor, 24" PU mag rear wheel



KJW-634E

- 1) Metal steel frame
- 2) Flip-up armrest and detachable elevating footrest
- 3) 8" PVC front castor, 24" PU mag rear wheel



KJW-653(ZE)

- 1) Metal steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 200*50 PVC front castor, 24" PU mag rear wheel

型号	尺寸 (mm)	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高
KJW-611	60/65	23	41/46	59	20	49	91	104	40	43	100	81X23X80	16	18		
KJW-621	65	28	46	59	20	50	87	106	42	37	100	93X28X87	17	19		
KJW-634E	67	29	46	60	20	50	91	102	40	41	100	78X29X91	16	18		

型号	尺寸 (mm)	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高
KJW-636M	67	29	46	60	20	50	91	102	40	41	100	78X29X91	16	18		
KJW-635M	71	28	56	61	20	52	93	107	45	42	125	81X24X93	18	20		
KJW-653ZE	73	30	56	60	20	52	93	107	45	42	125	81X26X93	18	20		



KJW-675(B)

- 1) Metal steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 8" PVC front castor and 24" PU mag rear wheel



KJW-675(A)

- 1) Metal steel frame
- 2) Flip-up armrest and detachable footrest
- 3) 8" PVC front castor and 24" solid rear wheel



KJW-671 (HE)

- 1) Chromed metal steel frame
- 2) Detachable armrest and elevating footrest
- 3) 8" PVC front castor, 24" solid rear wheel
- 4) Reclining high back



KJW-676(B)

- 1) Metal or chromed steel frame
- 2) Fixed armrest and footrest
- 3) 8" PVC front castor and 24" solid rear wheel



KJW-673(L)

- 1) Aluminum frame, with united brake
- 2) Detachable armrest and footrest
- 3) 6" solid castor and 12" PU back wheel.



KJW-672(HE)

- 1) Chromed metal steel frame
- 2) Detachable armrest and elevating footrest
- 3) 8" PVC front castor, 24" solid rear wheel
- 4) Reclining high back

型号	尺寸 (mm)													
	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高	座宽	
KJW-675B	67	29	46	60	20	54	91	105	43	36	100	78X34X91	20	22
KJW-675A	65	26	46	59	20	54	91	105	43	36	100	78X32X91	20	22
KJW-671HE	69	32	46	59	20	54	124	113	43	45	100	91X34X93	21	23

型号	尺寸 (mm)													
	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高	座宽	座深	座高	座宽	
KJW-676B	63	23	46	59	20	54	87	106	43	33	100	93X26X87	19	21
KJW-673L	63	23	46	59	20	54	87	106	43	33	100	93X26X87	19	21
KJW-672HE	65	31	46	59	20	53	127	118	43	44	100	90X31X95	25	27



KJW-R502

- Desinged to give maximum mobility for the user
- Product has many important features.e.g.Quick adjust handle height and easy
- Fold system Superior manouverability and seating comfort
- Fold into a compact unit for ease of transportation
- The rollator is more compact than other units when folded
- Low maintainance
- Suitable for indoor and outdoor use

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	folded width	castor size	unit weight	weight capacity
Size										
KJW-R502	min80cm max93cm	66.5cm	57.5cm	24.5cm	43.5cm	57cm	22cm	20cm	7.5kg	136.0kg



KJW-R503

- Heavy duty construction-with a safe working limit of 200kg
- Easy height adjustment snap button system
- Strong woven nylon sling backrest,conform to the shape of the user back
- Light touch lock brakes perfect for arthritis suffer
- Extra wide double stitched seat,breathable nylon cover
- 20cm flat ridged wheels minimises wheel wobble and increases grip
- Shopping bag has 75% more storage space
- New nylon bag prevents visibility of valuable items.e.g.wallet,hand bag,keys
- Foldable for ease of transportation
- Suitable for indoor and outdoor use

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	folded width	castor size	unit weight	weight capacity
Size										
KJW-R503	min83cm max95cm	77cm	67cm	34cm	41cm	51cm	49cm	20cm	9.4kg	200.0kg



KJW-R500

- Heavy duty construction-with a safe working limit of 130kg
- New"light touch" lock brake system
- New 20cm flat ridged wheels to minimise wheel wobble and increase grip
- Easy fold rounded backrest for comfort
- Can be folded without removing the storage bag
- New nylon bag prevents visibility of valuable items.e.g.wallet,handbag,keys
- For indoor and outdoor use

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	folded width	castor size	unit weight	weight capacity
Size										
KJW-R500	min81cm max94cm	69cm	58cm	33cm	36cm	56.5cm	32cm	20cm	7.45kg	130.0kg



KJW-R501

- Heavy duty construction-with a safe working limit of 130kg
- New"light touch" lock brake system
- New 20cm flat ridged wheels to minimise wheel wobble and increase grip
- Easy fold rounded backrest for comfort
- Can be folded without removing the storage bag
- New nylon bag prevents visibility of valuable items.e.g.wallet,handbag,keys
- For indoor and outdoor use

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	folded width	castor size	unit weight	weight capacity
Size										
KJW-R501	min81cm max93cm	68cm	61cm	33cm	36cm	49cm	26cm	15cm	6.6kg	130.0kg



KJW-265L

- New aluminium oval tube frame construction
- Safe working load of 180kg
- The frame is widened between the rear wheels for maximum walking space
- The gutter arm walker is foldable for transportation
- The oval tube provides a directional height adjustment function
- (Press button adjustment)
- The handles bars are multi adjustment for angle and fore/aft function
- Suitable for indoor and outdoor use

Item No.	overall height	overall length	overall width	armrest length	folded width	castor size	unit weight	weight capacity
Size								
KJW-265L	min103.5cm max123cm	86cm	63.5cm	min35cm max45cm	22.5cm	20cm	9.2kg	180.0kg



KJW-266L

- New aluminium oval tube frame construction
- Safe working load of 200kg
- The frame is widened between the rear wheels for maximum walking space
- Extra depth gutter arm supports for maximum comfort and safety
- The oval tube provides a directional height adjustment function
- (Press button adjustment)
- The handles bars are multi adjustment for angle and fore/aft function
- Suitable for indoor and outdoor use

Item No.	overall height	overall length	overall width	armrest length	castor size	unit weight	weight capacity
Size							
KJW-266L	min102cm max134.5cm	84cm	66cm	min28cm max40.5cm	10cm	10.05kg	200.0kg



KJW-263L

- New aluminium circular tube frame construction
- Safe working load of 160kg
- The frame is widened between the rear wheels for maximum walking space
- The Gutter arm walker is foldable for transportation
- The circular tube provides a directional height adjustment function
- (Press button adjustment)
- The handles bars are multi adjustment for angle and fore/aft function
- Suitable for indoor and outdoor use



KJW-264L

- Aluminium frame
- Available in modern colours
- Packed in retail packaging
- Accessories are available for modifications

Item No.	overall height	overall length	overall width	armrest length	folded width	castor size	unit weight	weight capacity
Size								
KJW-263L	min106.5cm max129cm	75cm	61cm	min42cm max54cm	32cm	10cm	10.2kg	160.0kg

Item No.	overall height	overall length	overall width	unit weight	weight capacity
Size					
KJW-264L	min110.5cm max125cm	52cm	57.5cm	4.0kg	100.0kg



KJW-2506-5"



KJW-250L



KJW-254L



KJW-254L-KD



KJW-263L



KJW-261L-5"

Item No.	overall height	overall length	overall width	unit weight
Size				
KJW-2506-5"	77.5-90cm	59cm	48.5cm	14.5kg
KJW-250L	76-88.5cm	54cm	45cm	11kg
KJW-254L	74.5-87cm	55cm	50cm	12.5kg

Item No.	overall height	overall length	overall width	unit weight
Size				
KJW-254L-KD	76-88.5cm	55cm	48.5cm	11.0kg
KJW-263L	74.5-87cm	55cm	50cm	11.2kg
KJW-261L-5"	73-85.5cm	48cm	55cm	6.5kg



KJW-261



KJW-255L



KJW-262



KJW-253L



KJW-256L



KJW-254L-5"

Item No.	overall height	overall length	overall width	unit weight
Size				
KJW-261	90-116cm	71cm	55cm	5.8kg
KJW-255L	73-85.5cm	55cm	43cm	5.0kg
KJW-262	90-116cm	71cm	55cm	5.8kg

Item No.	overall height	overall length	overall width	unit weight
Size				
KJW-253L	73-85.5cm	55cm	43cm	5.0kg
KJW-256L	73-85.5cm	48cm	55cm	4.3kg
KJW-254L-5"	81-91cm	59cm	53cm	7.0kg



KJW-3602L

- This transfer bench has a safe working load of 250kg
- Oval tubing provides a strong and robust frame
- The backrest is adaptable to fit on either side of the bench
- Sure grip suction tips are used to anchor the position, preventing it from slipping
- Suitable for hospitals, nursing homes and other institutional use



KJW-3603L

- Constructed from aluminium oval tubing
- New "S" bend side frame provides maximum space for the patient to fit between the armrests
- Offers a safe working load of 280kg
- Super light aluminium frame
- Fitted with slip resistant suction tips
- With plastic backrest
- Suitable for hospitals, nursing homes and other institutional use
- Low maintenance
- Heavy duty

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
Size									
KJW-3602L	min83.5cm max96cm	55cm	78cm	41.5cm	68cm	min49cm max61.5cm	26.5cm	6.8kg	250kg

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
Size									
KJW-3603L	min89.5cm max99cm	54cm	68.5cm	34.5cm	50.2cm	min45cm max55cm	27.5cm	7.6kg	280kg



KJW-3502L

- Constructed from aluminium oval tubing
- New design provides maximum space for the patient to fit between the arm rests
- This product offers a higher safe working load
- Ultra light aluminium frame
- Fitted with slip resistant suction tips
- Painted up to 110 microns. That is 3x thicker than normal coating to prevent scratches
- Suitable for hospitals, nursing homes and other institutional use
- Low maintenance

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
Size									
KJW-3502L	min66cm max78.5cm	46cm	57cm	33.5cm	40cm	min44.5cm max57.5cm	21cm	2.9kg	205kg



KJW-3503L

- Constructed from aluminium oval tubing
- New "S" bend design provides maximum space for the patient to fit between the arm rests
- The seat is made from a high density plastic to minimise flex
- This product offers a higher safe working load
- Super light aluminium frame
- Fitted with slip resistant suction tips
- Painted up to 110 microns. That is 3x thicker than normal coating to prevent scratches
- Suitable for hospitals, nursing homes and other institutional use
- Low maintenance

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
Size									
KJW-3503L	min76.8cm max89.2cm	52.2cm	55.2cm	42cm	44cm	min42cm max54.3cm	22cm	3.15kg	140kg



KJW-3604L

- The new seating provides full access to all the body parts of the patient
- Absolutely no need for the user to stand off the chair to have their body parts washed
- The frame is made from aluminium oval tubing
- Offer a safe working load of 250kg
- Soft polyurethane seat is waterproof
- With plastic backrest
- Suitable for hospitals, nursing homes and other institutional use
- Low maintenance



KJW-3504L

- Constructed from aluminium tubing
- New "S" bend side frame provides maximum space for the patient to fit between the armrests
- Seat and the lid are made from a high density plastic to minimise flex
- Offer a higher safe working load
- Ultra light aluminium frame
- Fitted with rubber suction tips
- Painted up to 110 microns. That is 3x thicker than normal coating to prevent scratches
- Suitable for hospitals, nursing homes and other institutional use
- Low maintenance
- Heavy duty

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
Size									
KJW-3604L	min94.5cm max103.5cm	56cm	68.5cm	45cm	48cm	min49cm max59cm	19.3cm	8.85kg	250.0kg

Item No.	overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
Size									
KJW-3504L	min62cm max74.5cm	53cm	56.5cm	42.3cm	35.3cm	min42cm max55.3cm	19cm	3.4kg	205.0kg



KJW-3501LB



KJW-3600LB



KJW-3600LA



KJW-3501LA



KJW-3601LA



KJW-3500L

SIZE MODEL						PCS CN
KJW-3501LB	34.5cm	34.5cm	35-50cm	35-50cm	43X43X42cm	4
KJW-3600LA	47cm	45cm	40-52.5cm	73-85.5cm	65X37X50cm	4
KJW-3601LA	72cm	45cm	43.5-44.5cm	79-89cm	77X22X45cm	2

SIZE MODEL						PCS CN
KJW-3600LB	48cm	43cm	41-53.5cm	71.5-84cm	57.5X33X49cm	4
KJW-3501LA	47cm	45cm	40-52.5cm	40-52.5cm	52X39X33cm	4
KJW-3500L	76.5cm	24.5cm	-	18cm	79X288.5X37cm	5



KJW-334L



KJW-335L



KJW-336L



KJW-337L



KJW-332



KJW-331

SIZE MODEL						PCS CN
KJW-334L	50cm	45cm	36.5-46.5cm	61-71cm	41X16X68cm	2
KJW-335L	48.5cm	57cm	43-53cm	92-101cm	56X18X79cm	1
KJW-336L	57cm	49.5cm	41-51cm	81-91cm	56X18X75cm	1

SIZE MODEL						PCS CN
KJW-337L	56cm	46cm	47-54.5cm	91-98.5cm	56X18X79cm	1
KJW-332	52cm	44cm	43.5-53.5cm	70-80cm	53X30X70cm	2
KJW-331	51.5cm	48cm	46.5-56.5cm	72.5-82.5cm	53X30X73cm	2



KJW-333L



KJW-328L



KJW-324








KJW-321








KJW-322



KJW-323

SIZE MODEL						PCS CN
KJW-333L	51.5cm	44.5cm	37-47cm	62.5-72.5cm	50X30X70cm	2
KJW-328L	58.5cm	58cm	41.5-51.5cm	66.5-76.5cm	90X31X64cm	2
KJW-324	42cm	33cm	46.5cm	46.5cm	60X57X41cm	5

SIZE MODEL						PCS CN
KJW-321	37.5cm	37cm	40.5cm	40.5cm	80X38X26cm	6
KJW-322	45cm	36cm	42cm	83.5cm	73X48X38cm	4
KJW-323	46cm	48.5cm	45.5cm	77.5cm	78X39X55cm	4



KJW-314S



KJW-313S



KJW-312



KJW-315S

SIZE MODEL						PCS CN
KJW-314S	91cm	44cm	49cm	95cm	55X54X55cm	1
KJW-312	91cm	44cm	49cm	95cm	55X54X55cm	1

SIZE MODEL						PCS CN
KJW-313S	91cm	44cm	49cm	95cm	55X54X55cm	1
KJW-315S	82cm	46cm	55cm	104cm	55X54X55cm	1



KJW-171L(L, M, S)






KJW-171L(MA)



KJW-171L(MB)



KJW-171L(MC)

SIZE MODEL	HEIGHT		PCS CN	SIZE MODEL	HEIGHT		PCS CN	SIZE MODEL	HEIGHT		PCS CN
KJW-171S	96-116cm	101X33X29cm	10	KJW-171L	136-156cm	141X33X29cm	10	KJW-171L(MB)	98-158cm	65x28.5x34cm	10
KJW-171M	116-136cm	121X33X29cm	10	KJW-171L(MA)	119-139cm	65X28.5X34cm	10	KJW-171L(MC)	120-140cm	65X28.5X34cm	10



KJW-134LA



KJW-134LB



KJW-135L



KJW-132LA



KJW-132LB



KJW-131LB



KJW-131LA




KJW-121LA




KJW-121LC

KJW-121LB



SIZE MODEL	HEIGHT		PCS CN
KJW-134LA	87-175.5cm	80X26X29cm	20
KJW-134LB	87-175.5cm	80X26X29cm	20
KJW-135L	87-175.5cm	80X26X29cm	20
KJW-132LA	94-124.5cm	93X28X30cm	20
KJW-132LB	101-127cm	80X26X29cm	20

SIZE MODEL	HEIGHT		PCS CN
KJW-131LA	96.5-111.5cm	100X30X30cm	20
KJW-131LB	96.5-111.5cm	100X30X30cm	20
KJW-121LB	84cm	86X32X32cm	10
KJW-121LA	84cm	86X32X32cm	10
KJW-121LC	84cm	86X32X32cm	10



KJW-155L



KJW-162B



KJW-161



KJW-154L



KJW-165



KJW-128L



KJW-123LB



KJW-123LA



KJW-153L



KJW-162A



KJW-163



KJW-164S



KJW-164



KJW-120A



KJW-120LB



KJW-139L

SIZE MODEL	HEIGHT		PCS CN
KJW-155L	72-95cm	81X38X34cm	20
KJW-162B	94.5-104.5cm	112X50X22cm	10
KJW-161	72-95cm	81X38X34cm	20
KJW-154L	69.5-87cm	81X40X36cm	10
KJW-165	73.5-96.5cm	81X38X34cm	15
KJW-128L	72-95.5cm	78.5X37.5X33cm	20
KJW-123LB	80-90cm	62X16X28cm	20
KJW-123LA	80-90cm	62X16X28cm	20
KJW-153L	94.5-104.5cm	112X50X22cm	10
KJW-162A	66-89cm	73X29X32cm	10
KJW-163	73-96cm	73X29X32cm	10
KJW-164S	62-92cm	81X38X34cm	20
KJW-164	73-96cm	81X38X34cm	10
KJW-120A	70.5-93.5cm	68X16X27cm	20
KJW-120LB	70.5-93.5cm	68X16X27cm	20
KJW-139L	119cm	37X27X20cm	50



A01



A02



A03



A04



A05



A06



A07



A08



A09



A10



A11



A12



A13



A14



A15



A16



A17



A18



B01



B02



B03



B04



B05



B06



B07



B08



B09



B10



B11



B12



B13



B14



B15



B16



B17



B18



B19



B20



B21



B22



C01



C02




C03



C04



D01





GUANGDONG KAREWAY MEDICAL TECHNOLOGY CO.,LTD.

Add:One of No.19 jinshi Avenue,Nanhai National Eco-industrial Demonstration Zone,Danzao Town,Nanhai District,Foshan City,Guangdong,China.

Mob:0086-18823130539

Tel: 0086-757-85408801

Fax:0086-757-85409990

Website:www.kareway-med.com

E-mail:kareway01@kareway-med.com

kareway02@kareway-med.com