



KAREWAY

广东康进威医疗科技有限公司

01

广东康进威医疗科技有限公司是一家专业致力于医疗器械的研发，生产，销售于一体的高新科技企业。本公司的产品以医用护理设备，康复设备为主，主要包括电动护理床，手动护理床，电动轮椅，手动轮椅，拐杖，助行器等。本公司自成立以来，汇聚了一批专业素质过硬，充满活力的工程技术人员。我们一直坚持“诚信，创新，服务”为理念，提供一流的产品，为客户创造价值。我们除了坚持独特的研发技术外，还特别坚持完善的产品售后服务体系。

本厂通过ISO9001:2008、ISO13485、CE 认证。我们致力于用创新和人性化的设计为客户提供完全，实用，便捷的产品，期待与您的合作。



Guangdong kareway medical technology Co., LTD is a professional medical equipment manufacturer with its own R&D, production and sales facility in Foshan China. Our company's products covers a range of electric hospital beds, manual hospital beds, electric wheelchairs, manual wheelchairs, walkers, walking sticks and other rehabilitation devices; we have a very professional R&D team with advanced manufacturing skills and good quality control system, after sale service is also our focus.

Our company has been certified by ISO9001:2008, ISO13485:2003. and our products have been approved by CE certificate. Our vision is to provide the world class products and services to attract the customers from all over the world. we look forward to good cooperation with you.

KAREWAY



KAREWAY

广东康进威医疗科技有限公司

KAREWAY




QUALIFICATION





KJW-S331(G)

THREE CRANKS MANUAL ORTHOPEDICS BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;

Related parameters:

- Overall dimension :
- L2220xW920xH480mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;

Standard configuration:

- 3sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails
- 4pcs urine hooks;
- 1pc S.S I.V pole;
- 1set orthopedics traction frame;


Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



KJW-S432(G)

FOUR CRANKS MANUAL ORTHOPEDICS BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension :
- L2220xW920xH460-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between:
- 460-720mm;

Standard configuration:

- 4sets manual crank systems;
- control-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails
- 4pcs urine hooks;
- 1pc S.S I.V pole;
- 1set orthopedics traction frame;

Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



GUANGDONG KAREWAY MEDICAL TECHNOLOGY CO.,LTD.

Add:One of No.19 jinshi Avenue,Nanhai National Eco-industrial Demonstration Zone,Danzao Town,Nanhai District,Foshan City,Guangdong,China.

Mob:0086-18823130539

Tel: 0086-757-85408801

Fax:0086-757-85409990

Website:www.kareway-med.com

E-mail:kareway01@kareway-med.com

kareway02@kareway-med.com