

AmcareMed for your and your patient's Safety

CONVENTIONAL FLOWMETER

TESTED FUNCTIONALITIES

AmcareMed flowmeters are instant flow measuring devices suitable for the dosage of medical gases. They are manufactured in version with in single or twin configuration in order to allow a double and independent gas supply using a single gas source. Flowmeters shall be supplied with gas specific probes and only attach to the corresponding medical gas terminal unit. Use in a clinical environment. The pressure dome and accuracy is tested which ensures the safety of both user and patient. Calibrated at 50 psi inlet pressure.

AVAILABLE IN US AND ISO COLORS



FH-B-O-15

FH-D-O-15

FH-D-A-15

FH-A-O-15

FL-DB-O

Working pressure	4.5bar
Setting ranges	1 to 15 L/min 1 to 10 L/min 0.5 to 15 L/min
Accuracy of flow rate	± 10 %
Screw thread for accessory	M 10 * 1 or 1/8 " adapter
Material	Flowmeter Aluminum or Brass Humidifier Bottle Polycarbonate plastic



RTM3M

FLOWMETER REGULATORS

- Gas Type : Medical Oxygen , Medical Air use
- Flow range 0-15 LPM
- Back pressure compensated flow meters
- 3000 psi maximum inlet pressure
- Easy to read flow rate window
- Complies to CGA standards



AM-O-GR-582



AM-O-GR-02



AM-O-GR-04



AM-O-GR-06



AM-O-GR-10

AVAILABLE IN ANY CHOICE OF INLET CONNECTION



CGA992 (BS)
Handtight



Japanese
Handtight



CGA540



CGA870



Rotary CGA870



CGA992 (BS)

CLICK STYLE

- Click style is available in 50 psi preset
- Available with a Nut and Nipple or Yoke Style
- Available in a variety of inlet adapter
- Available in a variety of colors (green is standard)
- Aluminum body with brass core to ensure safety
- Flowmeter available in 0-4 LPM, 0-8LPM, 0-15LPM

CS-R-870-121	CS-R-870-131	CS-R-540-121	CS-R-540-131	CS-R-870-135
CGA870 connection Outlet DISS O2	CGA870 connection Outlet Hose Barb	CGA540 connection Outlet DISS O2	CGA540 connection Outlet Hose Barb	CGA870 connection Outlet Hose Barb & DISS

ACCESSARY



Medical Gas Hose

- Available in custom length and configurations
- Available in a variety of gases
- Available in USA and ISO color



Twin Schrader Adaptor

- Twin outlet valve for attaching two hoses normally to one regulator
- Individual control valves for each outlet
- 100% air tightness test



TS-A-XX

Rail Clamp

- Rail clamps for mounting onto equipment or for mounting equipment on to rail clamps
- Anodised aluminium



RC-50

Standard Regulator & Self Sealing Valve

- Bull Nose & Pin Index Cylinder Fitting
- Brass Chrome Plated Body
- Self Sealing Valve
- Available oxygen, medical air



AM-R-SS-01



AM-R-SS-02

Monitor Bracket

- Pallet left and right angular: each 180 ° (with damping)
- Pallet pitch adjustment angle: unidirectional maximum 15°
- Bracket load: maximum 40kg
- Can be customized according to the monitor model



AM-M-60

Cylinder Trolley

Cylinder trolley used to solve the transportation problem of hospital oxygen cylinders. It can be adjusted to a variety of heights, the body painted a layer of black paint. It is easy to use and fully adjustable.



CY-C-A



CY-C-B