

INSPIRAL



STAINLESS STEEL
HOSPITAL FURNITURE
2019



SK60.30



SK60.10

Scrub Sinks

- Made of AISI 304 stainless steel designed to be used by the surgery team for hand and nail cleaning before entering the sterile zone
- Photocell display screen that shows the temperature and duration of water alternatively from each faucet
- Specially designed with round concave corners, easy to clean and antibacterial material considering the areas of use & user's requests and to be compatible with international hygienic norms
- Photocell and knee controlled water flow and knee controlled soap flow at each station
- Thermostatic mixing valve for cold and hot water

Article Codes	L (mm)	W (mm)	H (mm)	Explanation
SK60.10	800	635	1180	Single
SK60.20	1600	635	1180	Double
SK60.30	2200	635	1180	Triple



Water temperature and timer can be monitored on LCD display

SK60.20



Faucets

- Water and soap faucets are designed to keep hand sanitation
- High grade chrome plating for long life span
- Equipped with water saver

Sink

- Designed to prevent splash during washing
- Easy to clean, hygienic surfaces