


MANUAL HOSPITAL BED

A photograph of a family celebrating a patient's recovery in a hospital room. A young girl with blonde hair is sitting up in a hospital bed, wearing a light blue hospital gown. She is smiling and looking towards a young boy sitting on the edge of the bed. The boy is wearing a striped shirt and is also smiling. A man in a purple polo shirt and a woman in a pink shirt are sitting on the bed with them, all looking at a card or document. A teddy bear is sitting on the bed next to the girl. There are balloons and flowers in the background, suggesting a birthday or a special occasion. The scene is set in a bright, clean hospital room with a window in the background.



KJW-S332(SN)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: With food table

Related parameters:

- Overall dimension: L2170xW980xH450-885mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 450-885mm;

Standard configuration :

- 3 sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 4 pcs S. S side rails. ;
- 4pcs urine hooks;
- 1pc S.S I.V pole;

Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-S308(LN-B)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2160xW920xH480-760mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 480-760mm;

Standard configuration:

- 3 set manual crank systems;
- 4 pcs luxurious castors, 4 pcs legs without castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4 pcs urine hooks;
- 1 pcs S.S.IV pole;



Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;



KJW-S303(PZ-C)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170 × W990 × H430–720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430–720mm;

Standard configuration:

- 3 sets manual crank systems;
- Central-controlled brake system;
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair Luxurious PP siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;
- Sundries basket;



KJW-S301(LN-C)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension: L2170xW990xH430-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430-720mm;

Standard configuration :

- 3 sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;
- Sundries basket;



KJW-S302(LZ-C)

THREE CRANKS MANUAL CARE BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:


- Overall dimension: L2170xW990xH430-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between: 430-720mm;

Standard configuration:

- 3 sets manual crank systems;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;

Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;
- Sundries basket;



KJW-S231(LN)

TWO CRANKS MANUAL CARE BED

Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:



- Overall dimension:
L2100 × W970 × H480mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;

Standard configuration :

- 2 sets manual crank systems ;
- 4pcs luxurious castors ;
- 1 pair plastic head/foot-end boards ;
- 1 pair Aluminium siderails
- 4pcs urine hooks ;
- 1pc S.S. I.V pole ;

Options:

- Medical mattress ;
- Dining board ;
- Bedside cabinet ;
- Sundries basket ;



KJW-S331(LN)

THREE CRANKS MANUAL CARE BED

Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension:
L2100 × W970 × H360-680mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between:
360-680mm;

Standard configuration :

- 3 sets manual crank systems ;
- 4pcs luxurious castors ;
- 1 pair plastic head/foot-end boards ;
- 1 pair aluminium siderails
- 4pcs urine hooks ;
- 1pc S.S. I.V pole ;


Options:

- Medical mattress ;
- Dining board ; Bedside cabinet ;



KJW-S331(G)

THREE CRANKS MANUAL ORTHOPEDICS BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;

Related parameters:

- Overall dimension :
- L2220xW920xH480mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;

Standard configuration:

- 3sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails
- 4pcs urine hooks;
- 1pc S.S I.V pole;
- 1set orthopedics traction frame;


Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



KJW-S432(G)

FOUR CRANKS MANUAL ORTHOPEDICS BED



Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

Related parameters:

- Overall dimension :
- L2220xW920xH460-720mm;
- Back rest tilting: $80^{\circ} \pm 5^{\circ}$;
- Knee board tilting: $40^{\circ} \pm 5^{\circ}$;
- Hi-lo adjustment between:
- 460-720mm;

Standard configuration:

- 4sets manual crank systems;
- control-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails
- 4pcs urine hooks;
- 1pc S.S I.V pole;
- 1set orthopedics traction frame;

Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



GUANGDONG KAREWAY MEDICAL TECHNOLOGY CO.,LTD.

Add:One of No.19 jinshi Avenue,Nanhai National Eco-industrial Demonstration Zone,Danzao Town,Nanhai District,Foshan City,Guangdong,China.

Mob:0086-18823130539

Tel: 0086-757-85408801

Fax:0086-757-85409990

Website:www.kareway-med.com

E-mail:kareway01@kareway-med.com

kareway02@kareway-med.com