

VACUUM REGULATOR

VR-I-W Vacuum Regulator

VR-I-W

- Available Gauges : 0-160 mmhg, 0-300 mmhg, 0-760 mmhg
- Vacuum trap prevents overflow of waster fluid from entering the suction regulator
- Lightweight and durable
- Large control knob for easy suction adjustable
- Color coded flow ranges easy to read
- Dual scale mmhg and kpa



VR-I-XX Vacuum Regulator


VR-I-XX

- Available Gauges : 0-400 mmhg, 0-760 mmhg
- Vacuum trap prevents overflow of waster fluid from entering the suction regulator
- Lightweight and durable
- Large control knob for easy suction adjustable
- Color coded flow ranges easy to read
- Dual scale mmhg and kpa



NURSE CALL SYSTEM

- Two way calling between patients and nursing staff via the hospital' s network
- With talk back function
- The sound of talk back between host and extension is clear
- Convenience to change the number of extension



Layout drawing



Mobile Trolley



Venturi Type Vacuum Regulator



Slide Bracket



1L Suction Jar



2L Suction Jar



4L Suction Jar



Host machine



Nurse call extension



List of patient



Toilet room extension



Display screen