


## KJW-S402(PZ)

### FOUR CRANKS MANUAL CARE BED



#### Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg;

#### Related parameters:

- Overall dimension: L2170 × W990 × H460–750mm;
- Back rest tilting:  $80^{\circ} \pm 5^{\circ}$  ;
- Knee board tilting:  $40^{\circ} \pm 5^{\circ}$  ;
- Hi-lo adjustment between: 460–750mm;
- Reverse and trendelenburg  $15^{\circ} \pm 2^{\circ}$  ;

#### Standard configuration :

- 4 sets manual crank systems;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


#### Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;



## KJW-S401(LN-C)

### FOUR CRANKS MANUAL CARE BED



#### Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;
- D: Trendelenburg and Reverse Trendelenburg;

#### Related parameters:

- Overall dimension: L2170 × W990 × H430–720mm;
- Back rest tilting:  $80^{\circ} \pm 5^{\circ}$  ;
- Knee board tilting:  $40^{\circ} \pm 5^{\circ}$  ;
- Hi-lo adjustment between: 430–720mm;
- Reverse and trendelenburg  $15^{\circ} \pm 2^{\circ}$  ;

#### Standard configuration:

- 4 sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;


#### Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;



## KJW-S331(G)

### THREE CRANKS MANUAL ORTHOPEDICS BED



#### Main functions:

- A: Back section up/down;
- B: Knee section up/down;

#### Related parameters:

- Overall dimension :
- L2220xW920xH480mm;
- Back rest tilting:  $80^{\circ} \pm 5^{\circ}$  ;
- Knee board tilting:  $40^{\circ} \pm 5^{\circ}$  ;

#### Standard configuration:

- 3sets manual crank systems;
- 4pcs luxurious castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails
- 4pcs urine hooks;
- 1pc S.S I.V pole;
- 1set orthopedics traction frame;


#### Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



## KJW-S432(G)

### FOUR CRANKS MANUAL ORTHOPEDICS BED



#### Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

#### Related parameters:

- Overall dimension :
- L2220xW920xH460-720mm;
- Back rest tilting:  $80^{\circ} \pm 5^{\circ}$  ;
- Knee board tilting:  $40^{\circ} \pm 5^{\circ}$  ;
- Hi-lo adjustment between:
- 460-720mm;

#### Standard configuration:

- 4sets manual crank systems;
- control-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair aluminium siderails
- 4pcs urine hooks;
- 1pc S.S I.V pole;
- 1set orthopedics traction frame;

#### Options:

- Medical mattress;
- Dining board; Bedside cabinet;
- Sundries basket;



## **GUANGDONG KAREWAY MEDICAL TECHNOLOGY CO.,LTD.**

Add:One of No.19 jinshi Avenue,Nanhai National Eco-industrial Demonstration Zone,Danzao Town,Nanhai District,Foshan City,Guangdong,China.

Mob:0086-18823130539

Tel: 0086-757-85408801

Fax:0086-757-85409990

Website:[www.kareway-med.com](http://www.kareway-med.com)

E-mail:[kareway01@kareway-med.com](mailto:kareway01@kareway-med.com)

[kareway02@kareway-med.com](mailto:kareway02@kareway-med.com)