




## KJW-S202(PN-C)

### TWO CRANKS MANUAL CARE BED



#### Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

#### Related parameters:

- Overall dimension: L2170 × W990 × H480mm;
- Back rest tilting:  $80^{\circ} \pm 5^{\circ}$  ;
- Knee board tilting:  $40^{\circ} \pm 5^{\circ}$  ;

#### Standard configuration :

- 2 sets manual crank systems;
- Central-controlled brake system,
- 4pcs luxurious central locking castors;
- 1 pair plastic head/foot-end boards;
- 1 pair L type siderails;
- 4pcs urine hooks;
- 1pc S.S I.V pole;



#### Options:

- Medical mattress;
- Dining board;
- Bedside cabinet;
- Sundries basket;



## KJW-S201(LN-C)

### TWO CRANKS MANUAL CARE BED



#### Main functions:

A:Back section up/down;  
B:Knee section up/down;

#### Related parameters:


Overall dimension : L2170xW990xH480mm;  
Back rest tilting:  $80^{\circ} \pm 5^{\circ}$  ;  
Knee board tilting:  $40^{\circ} \pm 5^{\circ}$  ;

#### Standard configuration :

2 sets manual crank systems ;  
4pcs luxurious castors ;  
1 pair plastic head/foot-end boards ;  
1 pair L type siderails ;  
4pcs urine hooks ;  
1pc S.S I.V pole ;

#### Options:

Medical mattress ;  
Dining board ;  
Bedside cabinet ;  
Sundries basket ;



## KJW-S231(LN)

TWO CRANKS MANUAL CARE BED

### Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

### Related parameters:



- Overall dimension:  
L2100 × W970 × H480mm;
- Back rest tilting:  $80^{\circ} \pm 5^{\circ}$  ;
- Knee board tilting:  $40^{\circ} \pm 5^{\circ}$  ;

### Standard configuration :

- 2 sets manual crank systems ;
- 4pcs luxurious castors ;
- 1 pair plastic head/foot-end boards ;
- 1 pair Aluminium siderails
- 4pcs urine hooks ;
- 1pc S.S. I.V pole ;

### Options:

- Medical mattress ;
- Dining board ;
- Bedside cabinet ;
- Sundries basket ;



## KJW-S331(LN)

THREE CRANKS MANUAL CARE BED

### Main functions:

- A: Back section up/down;
- B: Knee section up/down;
- C: Whole bed up/down;

### Related parameters:

- Overall dimension:  
L2100 × W970 × H360-680mm;
- Back rest tilting:  $80^{\circ} \pm 5^{\circ}$  ;
- Knee board tilting:  $40^{\circ} \pm 5^{\circ}$  ;
- Hi-lo adjustment between:  
360-680mm;

### Standard configuration :

- 3 sets manual crank systems ;
- 4pcs luxurious castors ;
- 1 pair plastic head/foot-end boards ;
- 1 pair aluminium siderails
- 4pcs urine hooks ;
- 1pc S.S. I.V pole ;

### Options:

- Medical mattress ;
- Dining board ; Bedside cabinet ;



## **GUANGDONG KAREWAY MEDICAL TECHNOLOGY CO.,LTD.**

Add:One of No.19 jinshi Avenue,Nanhai National Eco-industrial Demonstration Zone,Danzao Town,Nanhai District,Foshan City,Guangdong,China.

Mob:0086-18823130539

Tel: 0086-757-85408801

Fax:0086-757-85409990

Website:[www.kareway-med.com](http://www.kareway-med.com)

E-mail:[kareway01@kareway-med.com](mailto:kareway01@kareway-med.com)

[kareway02@kareway-med.com](mailto:kareway02@kareway-med.com)