

Electric obstetric delivery bed series

KJW-OP-D9400

Electric Obstetric table

KJW-OP-D9400 electric delivery bed up-down, back up-down, seat section up-down (or Trendelenburg and reverse-trendelenburg) are all controlled by

- The mattress is seamless foaming pad, which is beautiful and easy to clean.
- The leg board can be opened 90°, which is convenient for operation.
- The foot controller is standard configuration.



Main specification

Length/Width	1900mm/600mm	Back Section Up/Down	≥70°/≥10°
Height	600mm-900mm	Leg section opened	≥90°/≥10°
Trendelenburg/Revers Trendelenburg	≥20°/≥20°	Load Capacity	250kg



All the cover, able top, basin are all made of stainless steel 304#, which is convenient for clean and disinfection.

The examination light is optional, which can be used for gynecology examination.

Hydraulic obstetric delivery bed

KJW-OP-S8400
Hydraulic Obstetric table

KJW-OP-S8400 delivery bed table up-down is hydraulic controlled.

The column and base cover, basin are all made of stainless steel 304#, which is convenient for clean and disinfection.



Main specification

Length/Width	2000mm/600mm
Height	620mm/860mm
Trendelenburg/Revers Trendelenburg	≥20°/≥20°
Back Section Up/Down	≥65°/≥10°
Leg Section Down	250kg

Manual obstetric delivery bed

KJW-OP-S8200
Hydraulic Obstetric table



Main specification

Length/Width/Height	1900mm/600mm/750mm
Back Section Up/Down	≥45°/≥10°
Seat Section Up/Down	≥25°/≥8°
Leg Section Up/Down	≥25°/≥25°

KJW-OP-S8200 Ordinary obstetric bed is used for obstetrics gynecology examination. The bed top and main structure are all made of stainless steel 304#



GUANGDONG KAREWAY MEDICAL TECHNOLOGY CO.,LTD.

Add:One of No.19 jinshi Avenue,Nanhai National Eco-industrial Demonstration Zone,Danzao Town,Nanhai District,Foshan City,Guangdong,China.

Mob:0086-18823130539

Tel: 0086-757-85408801

Fax:0086-757-85409990

Website:www.kareway-med.com

E-mail:kareway01@kareway-med.com

kareway02@kareway-med.com