




## 型号: KJW-D301(PZ-C)

品名: 三功能电动护理床

THREE FUNCTIONS CARE BED



KJW-D301(LZ-C)

### 主要功能:

- A:背部升降
- B:腿部升降
- C:整体升降

### 产品参数:

床体规格: L2170 × W1030 × H460-720mm  
 背部倾斜:  $80^{\circ} \pm 5^{\circ}$   
 腿部倾斜:  $40^{\circ} \pm 5^{\circ}$   
 床体可调高度: 460-720mm

### 标准配置:

高级医用电机三个, 控制盒一个;  
 控制手柄一个, 中控刹车系统;  
 豪华中控轮四个, 塑料床头尾板一对;  
 PP护栏一对, 引流挂钩四个;  
 不锈钢输液杆一支;

### 可选配:

后配电源UPS; 快速放倒功能CPR;  
 医用床垫; 餐板; 床头柜; 病历卡;  
 杂物篓。

### Main functions:

- A:Back section up/down.
- B:Knee section up/down.
- C:Whole bed up/down.

### Related parameters:

Overall dimension: L2170 × W1030 × H460-720mm  
 Back rest tilting:  $80^{\circ} \pm 5^{\circ}$   
 Knee board tilting:  $40^{\circ} \pm 5^{\circ}$   
 Hi-lo adjustment between: 460-720mm

### Standard configuration:

3pcs Advanced medical actuators;  
 1pcs control box, 1pcs handset;  
 central-controlled brake system;  
 4pcs luxurious central locking castors;  
 1 pair plastic head/foot-end board;  
 1 pair PP siderails;  
 4pcs urine hooks; 1pcs S.S. I.V pole;

### Options:

Backup power system (UPS);  
 CPRfunction; Medical mattress;  
 Diningboard; Bedside cabinet;  
 Patient file panel; Sundries basket.