

## INFINITUM SERIES

### IN 42 ELECTRONIC INTENSIVE CARE BED

INFINITUM hospital beds are a superior quality hospital bed that suits various environments - from med-surg to high acuity care. It is designed to maximize infection control and fall prevention, while being easy to use. With wide range of accessory options INFINITUM covers comprehensive needs of users.



- Electrical operated backrest, height and legrest adjustment
- Electrical operated Trendelenburg, Reverse Trendelenburg, fowler and vascular position
- Patient handset
- Foot end nurse control unit
- Footrest adjustment by ratchet
- Auto-CPR
- Auto Regression system on backrest
- Auto contour
- Hygienic design (PP), lockable and tuck away side rails
- Easy removable ABS mattress platform
- Removable head and footboards
- Electrostatic epoxy painted metal frame
- Height adjustable stainless steel IV pole
- Plastic protective crash bumpers on four corners
- Central Lockable Castors, 150mm
- Angle gauge for back
- Angle gauge for Trendelenburg and Reverse Trendelenburg positions
- Cardiac chair position with single button
- Shock position with single button
- Semi-Fowler position with single button
- Bed exit position with single button
- Examination position with single button
- Electronic CPR with one touch button
- Activation key at remote control to prevent any undesired positions
- Urine bag holder
- Urine basket (optional)
- Dual side manual CPR levers at backrest
- Rechargeable Battery Backup
- Bed size extension (optional)
- X-ray cassette holder at backrest (optional)
- IPX6 standard electronic svstem

www.optium.com.tr

### TECHNICAL DETAILS

Overall Length:	216 cm
Overall Width:	102 cm
Height:	37-80 cm
Trendelenburg:	0°-16°
Reverse Trendelenburg:	0°-16°
Backrest Angle:	0°-70°
Legrest Angle:	0°-40°
Castor Diameter:	15 cm
Under bed Clearance:	15 cm
Product Weight:	115 kg
Safe Working Load:	250 kg

### ELECTRICAL DATA

Power In:	1.5A max. At 230V AC 50-60 Hz
Electrical Shock Protection:	Class I, Type B
Liquid Ingress Protection:	IPX6
Batteries:	2x 12V Series Connected, sealed, rechargeable, Lead / Acid gel, 1.2 amp hrs.




**Accessories**



Central Brake Castors  
150mm



Foot-end Nurse  
Control Panel



Rechargeable Battery  
Backup



Dual Sided Manual  
CPR at Backrest



X-Ray Cassette Holder  
at Backrest



X-Ray Translucent HPL  
Backrest



Stationary Headboard



Control Panels Embedded  
on Side Rails



Urine Bag Basket



Bed Extension



Lifting Pole




Monitor & Defibrillator  
Shelf



Oxygen Cylinder Holder



Nurse Panel Shelf



Orthopedic Traction  
System

**Color Options**



P01



P02



P03



P04



P05



P06



P07



P08

www.optium.com.tr

